KULLANIM KILAVUZU

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 Shenzhen Homed Medical Device Co., Ltd. 3rd Floor, Block 2 Longgu Industrial Zone, Dalang Street, District, Longhua New District, Shenzhen 518109, People's Republic of China
www.homedgroup.com

 Shanghai International Holding Corp. GmbH

 0413

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lifecek

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USERS MANUAL

Please read the instructions carefully before operating the nebulizer.

1- What is the nebulizer?

The Lifechek JLN2308AS Plus Compressor Nebulizer System is a medical device designed to deliver medication prescribed by a physician to the airways of the lungs. The nebulizer and air compressor combined provides an aerosol mist for inhalation by the patient for medical treatment only.

2- How does the nebulizer work?

Lifechek JLN2308AS Plus Compressor Nebulizer System is compact, simple to operate and convenient to carry. It uses a non-oil lubricating, single-cylinder piston pump to create air pressure.

Air flows through the air filter into the Nebulizer's air compressor unit, out the jet outlet and then through the Nebulizer.

The main airflow into and out of the air compressor unit is illustrated in Figure 1 below.


Figure 1 - Working Principle


3- Contents

Before using this product, make sure you have the following items:

- 1 Pc. Main device (compressor)
- 1 Pc. Aerosol nebulizer
- 1 Pc. Aerosol tubing
- 1 Pc. Mask (adult & child)
- 1 Pc. Mouth piece

4- Components

The Standard Nebulizer components are diagramed below in.

Figures 2

5- How do I insert medication into the nebulizer?

To insert medication into the nebulizer, do the following.


5¹ Open the cap. (see Figure 3)

5² Fill the cup with medication as prescribed by your physician. Be sure that the medication does not exceed the liquid site line of 6 ml (6 cc).

5³ Close the cap.

5⁴ Connect the mouthpiece.

NOTE: The mouthpiece is to be used for oral inhalation only.

Figures 3

6- How do I operate the nebulizer?

To operate the nebulizer, do the following:

- 6¹** Plug in the air compressor but do not connect the nebulizer. 2- Make sure that the air filter is attached to the air compressor. 3- Turn on the power switch.
- 6²** After the air compressor runs properly for 1-2 minutes, connect the jet outlet of the air compressor to the inlet of the nebulizer by the air hose.
- 6³** Begin treatment according to your physician's instructions.
- 6⁴** Inhale and exhale naturally. Be sure to stop when the medication is completed.
- 6⁵** When you have completed use, turn the air compressor unit off before you disconnect the air hose.

WARNING: Do not use near flammable objects or locations near explosive materials.

7- What do I do after I've used the nebulizer?

To keep your nebulizer working properly, after each use do the following:

- 7¹** After you have disconnected the air hose, disconnect the mouthpiece.
- 7²** Open the cap and empty the cup. Wash the nebulizer with running water or soak in warm water for 15 minutes. For a more thorough, hygienic cleaning, add some vinegar in the water.

NOTE: Do not boil the nebulizer. Doing so could damage it.

- 7³** Dry the nebulizer completely before storing it.

8- How do I replace the air filter?

Generally, the air filter should be replaced every three months or after 300 uses. To replace the airfilter:

8¹ Remove the filter cover.

8² Replace the filter.

8³ Reattach the filter cover.

NOTE: Whenever the filter becomes wet, allow it to dry completely before using.

9- What maintenance for storage does the nebulizer require?

- Always handle your nebulizer with care.

- Disconnect the power plug from the wall outlet when not in use.

- Keep your nebulizer out of children's reach.

- If you store your nebulizer, do not bend the air hose.

- If you store your nebulizer, try to keep it in the following environmental ranges:

Temperature: -40 °C - 55 °C

Relative humidity: < 95%

Atmospheric pressure: 0.05 MPa 0.106MPa

-If possible, store your nebulizer in a well ventilated room.

10- Specification

Power: 220 VAC, 50 Hz

Power consumption: Below 65W

Operation time: 30 mins on/30 mins off

Nebulizer flow: ≥ 4 L

Compressor flow: ≥ 8 L

Max. medication capacity: 8 mL

Particle size: 0.5 - 6 µm

MMAD: 2.44 µm

Average mobilization rate: ≥ 0.2 ml/min.

Remain amount: < 1.2 ml

Noise level: Below 60dBA

Compressor pressure range: 21.75 - 58 Psi/150 - 400 Kpa/1.5 - 4.0 Bar

Operation pressure range: 11.6 - 21.7 Psi/80 - 150 Kpa/0.8 - 1.5 Bar

Operating temperature range: 10 °C - 50 °C (50 °F - 104 °F)

Operation humidity: 10% - 95% RH

Storage temperature range: -20 °C - 70 °C (-4 °F - 158 °F)

Storage humidity range: 10% - 95% RH

Dimension (L x W x H): 205 x 155 x 113 mm

Products weight: 1.55 Kg

Standard accessories: Compressor nebulizer machine

Adult mask

Child mask

Filter

Durable air tube

Medical cup (8cc)

Mouth piece

11- Periodic maintenance

The product does not require any periodic maintenance.

Information about usage

- All connections must be made carefully.
- The nebulizer should not be exposed to impacts.
- Usermanual and usage instruction should be read and applied.
- The nebulizer should not come into contact with any liquid while being used.
- Damage to the device may occur if it is not used under proper usage conditions.

12- Guidance and manufacturer's declaration

The JLN-2308AS is intended for use in the electromagnetic environment specified below. The customer or the user of the JLN-2308AS should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The JLN-2308AS uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The JLN-2308AS is intended for use in the electromagnetic environment specified below. The customer or the user of the JLN-2308AS should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<i>Electrostatic discharge (ESD) IEC 61000 4-2</i>	<i>±6 kV contact ±8 kV air</i>	<i>±6kV Contact ±8kV Air</i>	<i>Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. If ESD interfere with the operation of equipment, counter measurements such as wrist strap, grounding shall be considered.</i>
<i>Electrical fast transient/burst IEC 61000-4-4</i>	<i>±2 kV for power supply lines ±1 kV for input/output lines</i>	<i>±2 kV for Power supply lines</i>	<i>Mains power quality should be that of a typical commercial or hospital environment.</i>
<i>Surge IEC 61000 -4-5</i>	<i>±1 kV differential mode ±2 kV common mode</i>	<i>±1kV differential mode ±2kV common mode</i>	<i>Mains power quality should be that of a typical commercial or hospital environment.</i>
<i>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</i>	<i><5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec</i>	<i><5% UT for 0.5 cycle 40% UT for 5 cycles 70% UT for 25 cycles <5% UT for 5 s</i>	<i>Mains power quality should be that of a typical commercial or hospital environment. If the user of the JLN-2308AS requires continued operation during power mains interruptions, it is recommended that the JLN-2308AS be powered from an uninterruptible power supply or a battery.</i>
<i>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</i>	<i>3 A/m</i>	<i>3 A/m</i>	<i>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</i>

Guidance and manufacturer's declaration – electromagnetic immunity

The JLN-2308AS is intended for use in the electromagnetic environment specified below. The customer or the user of the JLN-2308AS should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the JLN-2308AS, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2$ 80 MHz to 800 MHz $d = 2.3$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3V/	

Recommended separation distances between Portable and mobile RF communications equipment and the JLN-2308AS

The JLN-2308AS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the JLN-2308AS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the JLN-2308AS as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2$	80 MHz to 800 MHz $d = 1.2$	800 MHz to 2.5 GHz $d = 2.3$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

14- Troubleshooting guide

Symptom	Probable Cause	Solution
Too much noise	No filter in the unit	- Put a filter in the unit.
Nebulisation may be distorted or there is no nebulisation	Main compressor switch is off.	- Turn the unit on.
	Nebulizer was not cleaned after its previous use.	- Clean the unit.
	Air hose is bent.	- Remove any bends or kinks in the air hose.
	The filter is clogged.	- Replace the filter.
	No medication is left.	- Add the appropriate amount of medication prescribed by your physician to the cup.
	The cap is not attached properly.	- Reattach the cap.
Water Drops have formed on the air hose.	The medication was over filled.	- Add appropriate amount of medication to cup, connect the air hose to the air compressor, and turn the unit on.
	The nebulizer was washed but not dried	- Open and close the nebulizer repeatedly with your finger until the drops disappear.

15- Information about usage

- All connections must be made carefully.
- The nebulizer should not be exposed to impacts.
- Usermanual and usage instruction should be read and applied.
- The nebulizer should not come into contact with any liquid while being used.
- Damage to the device may occur if it is not used under proper usage conditions.

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Shenzhen Homed Medical Device Co., Ltd. 3rd Floor, Block 2 Longgu Industrial Zone, Dalang Street, District, Longhua New District, Shenzhen 518109, People's Republic of China
www.homedgroup.com

EC REP Shanghai International Holding Corp. GmbH

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lifecek
is a **BEK** brand.